



Quantum

Cannulated Compression Screw

INTRODUCTION

Overview

Quantum Cannulated Compression Screw system offers the two kinds(2.2, 3.0) of screw head size to treat variety indication.

Continuous discussion with surgeons and personnel working in the OR and solution-finding efforts thereafter have enabled us to develop the Quantum Cannulated Compression Screw.

With this innovative product, we have been able to fulfill all of our valued customers' demands, ensuring the maximum quality and functionality level of the products we have to offer.

INDICATION

Fixation of fracture, osteotomies, fusions and arthrodeses



INSTRUMENT

I-Handle

D8P.01010

Screw Driver Shaft

Drill Guide

D8P.02010

K-Wire

K-Wire Depth Gauge

Drill Bit



For 2.2



D8P.22010

For 3.0



D8P.22020

For 2.2 / 3.0



For 0.8 (Φ2.2 SCREW)



D8P.04080

For 1.1 (Φ3.0 SCREW)



D8P.04110

For 0.8 K-Wire



D8P.03090

For 1.1 K-Wire



D8P.03120

For 2.2



D8P.31170


For 3.0



D8P.31210

Screws system

- ▶ Two head sizes accommodate patient anatomical condition.
- ▶ Wide range of screw lengths covers variety indication.
- ▶ High precise screw thread from accurate manufacturing process.
- ▶ Screws are made of most biocompatible material of titanium.
- ▶ Star shaped head design secures stable insertion with enough strength.



Length	Ø2.2				Ø3.0			
	Thread Length	Part No.	Thread Length	Part No.	Thread Length	Part No.	Thread Length	Part No.
	Headless Short		Headless Long		Headless Short		Headless Long	
6mm								
7mm								
8mm								
9mm								
10mm	4	D1D.22100			4	D2D.30100		
11mm	5	D1D.22110			4	D2D.30110		
12mm	5	D1D.22122			4	D2D.30120		
13mm	5	D1D.22130			4	D2D.30130		
14mm	5	D1D.22140			4	D2D.30140		
15mm	5	D1D.22150			4	D2D.30150		
16mm	5	D1D.22160			4	D2D.30160	5	D2E.30160
17mm	5	D1D.22170			4	D2D.30170	6	D2E.30170
18mm	5	D1D.22180			4	D2D.30180	6	D2E.30180
19mm	5	D1D.22190			4	D2D.30190	7	D2E.30190
20mm	5	D1D.22200			4	D2D.30200	7	D2E.30200
21mm	5	D1D.22210			4	D2D.30210	8	D2E.30210
22mm	5	D1D.22220	8	D1E.22220	4	D2D.30220	8	D2E.30220
23mm	5	D1D.22230			4	D2D.30230	8	D2E.30230
24mm	6	D1D.22240	8	D1E.22240	5	D2D.30240	8	D2E.30240
25mm	6	D1D.22250			5	D2D.30250	8	D2E.30250
26mm	6	D1D.22260	8	D1E.22260	5	D2D.30260	10	D2E.30260
27mm	6	D1D.22270			6	D2D.30270	10	D2E.30270
28mm	6	D1D.22280	9	D1E.22280	6	D2D.30280	10	D2E.30280
29mm	6	D1D.22290			6	D2D.30290	10	D2E.30290
30mm	6	D1D.22300	10	D1E.22300	7	D2D.30300	12	D2E.30300
32mm			11	D1E.22320	7	D2D.30320	12	D2E.30320
34mm			12	D1E.22340	8	D2D.30340	14	D2E.30340
35mm								
36mm			13	D1E.22360	9	D2D.30360	14	D2E.30360
38mm			14	D1E.22380	9	D2D.30380	16	D2E.30380
40mm			15	D1E.22400	10	D2D.30400	16	D2E.30400
42mm								
44mm								
45mm					15	D2D.30450	21	D2E.30450
46mm								
48mm								
50mm					20	D2D.30500	26	D2E.30500



OSTEONIC Co., Ltd.

1206Ho, 38, Digital-ro 29-gil, Guro-gu, Seoul, Korea P.O. #152-779
TEL: +822 6902 8400 FAX: +822 6902 8401 www.osteonic.com

CS - E - 1505